

## CYP0048- 2xTaq PCR MasterMix (with dye)

Cat#	Size
CYP0048	10*1ml

**Storage: 2-8°C**

**Attention:**

•This product contains phenol red, which may change color from orange to purple. **Color change does not affect product quality or use.**

**Product Description:**

- This product contains Taq DNA polymerase, dNTPs, MgCl<sub>2</sub>, and reaction buffer at a concentration of 2×. It has the advantages of being fast, easy, sensitive, specific, stable, and can minimize human errors. Only DNA template and primers are needed for use.
- This product is convenient and easy to use, and can avoid contamination during PCR operation. It only needs to add an appropriate amount of 2×Taq PCR MasterMix solution, add the template and primers, and add deionized water to make the MasterMix solution concentration 1× for the reaction. Please ensure that it is fully dissolved and mixed before use, and the operation should be carried out on ice.

**Product Contents:**

•Taq DNA polymerase (recombinant): 0.1 units/μl; MgCl<sub>2</sub>: 4 mM; dNTPs (dATP, dCTP, dGTP, dTTP): 0.4 mM

**Quality Control:**

•No exogenous nuclease activity was detected; it can effectively amplify single-copy genes in the human gene; no significant activity change was observed after one week of storage at room temperature.

**Scope of Application:**

- Gene detection: This product has small differences between batches, especially suitable for large-scale gene detection, semi-quantitative PCR experiments, and the detection of trace DNA.
- Amplification of DNA and some complex templates with special structures and high GC content (>60%) and secondary structures: amplification of DNA fragments, DNA labeling, primer extension, and sequence determination. PCR products are A-tailed and can be directly cloned into T/A vectors after purification.

**Reaction Examples: (Note: The following examples are only for reference. Actual reaction conditions may vary depending on the structure of the template, primers, etc. Adjust the reaction conditions accordingly.)**

•Use 2×Taq PCR MasterMix product to amplify a 1 kb fragment as a template with human genome DNA, and the reaction system is 25 μl. (If the reaction system is different, it can be increased or decreased in proportion.)

Template	10 pg-1 μg
Primer 1	0.1-1 μl
Primer 2	0.1-1 μl
2×MasterMix	12.5 μl
ddH <sub>2</sub> O	make up to 25 μl

**•PCR reaction cycle setting:**

94°C: 3 min	} 30 cycles
94°C: 30 sec	
55°C: 30 sec	
72°C: 1-2 kb/min	
72°C: 5 min	

**•Result detection:** After the reaction is finished, take 5 μl of the reaction product for agarose gel electrophoresis. If the product is amplified with a dye mix, it can be directly loaded onto the gel without adding sample buffer.

**This product is furnished for LABORATORY RESEARCH USE ONLY.  
 Not for diagnostic or therapeutic use.**